

Multiplex made effortless

the technology behind **Mabtech EYRA™**



Berend Bassleer
Product Manager, Instruments



Tyler Sandberg, PhD
Product Manager, Immunoassays




Pol del Aguila Pla
AI and Image Scientist, Qamcom

Our assays

We develop optimal monoclonal antibody pairs for sandwich assays

ELISpot


IFN- γ



Each spot is the footprint
of a responding cell

FluoroSpot


IL-22, IFN- γ , IL-5, IL-17A



Multiplexed
ELISpot

FociSpot


Influenza foci



Virus-
infected cells

ELISA

IFN- γ



Analyte
concentration

EYRAplex

29+ analytes



Bead-based
multiplex




IRIS 2 | ELISpot, FluoroSpot, FociSpot
ASTOR 2 | ELISpot




The Mabtech multiplex solution

A new bead-based multiplex platform – assay and instrument



EYRAplex




Mabtech EYRA™





The assay – EYRAplex

Magnetic beads with optimized antibody pairs and protein standards




The assay – EYRAplex

30+ distinct bead IDs and corresponding PE signals




EYRAplex step-by-step


Anyone can multiplex




1. Add capture bead mix




2. Add samples + standards



3. Add Detection mAb mix



4. Add SA-PE



5. Data acquisition

The assay – EYRAplex

The mAb pairs you trust, now in multiplex

Sample type and volume:
25 µl
serum, plasma, cell supernatant





- Kit includes:
 - Standard mixes
 - Capture bead mixes
 - Detection mAb mixes
 - Streptavidin-PE
 - Assay buffer, diluents, wash buffers
 - Black 96-well plate with lid and adhesive covers

Designed for **EYRA** and works
with most flow cytometers




The EYRAplex panels

A broad range of inflammatory and immune monitoring panels




- Th1/Th2 panels: 7-, 10-, and 18-plex
- Th17/Th22 panel: 13-plex
- Proinflammatory panels: 8- and 12-plex
- T reg panel: 12-plex
- *In vitro* cytokine release panel: 8- and 15-plex
- Chemokine panel: 6-plex
- Cytokine panel: 29-plex

- Cytokine storm: 20-plex
- Viral panel: 13- and 18-plex



How does EYRAplex measure up?

Highly sensitive mAb pairs make highly sensitive assays











EYRAplex on par and often more sensitive than competitors



EYRAplex assays are validated on EYRA and flow cytometry

Sensitivity and accuracy no matter what

- For the smoothest and fastest data acquisition, we recommend **EYRA**
 - **Insert plate, press read, and enjoy your data**
- EYRAplex kits perform just as well with most flow cytometers
 - Red and yellow-green/blue lasers needed



The instrument – Mabtech EYRA


Multiplex made effortless



The instrument – Mabtech EYRA

New, simple analysis of bead-based multiplex assays


- **Confocal microscopy**-based acquisition
- Bead complexes are imaged from below after excitation with the laser system
- RAWsphere™ algorithm ensures highly accurate bead identification and corresponding PE signal



Mabtech EYRA: a fluidics-free experience


No headaches related to the maintenance of a flow-based system

- Reading EYRAplex on EYRA means:
 - Sample stays in the well; **no-contact analysis**
 - **Worry-free experience**, no babysitting
 - No clogging, no sheath fluid or waste
 - **Plug-and-play** experience
 - True walk away acquisition, **~15 min read time**
 - **No start-up or shutdown procedures**



Growing plex numbers with Mabtech EYRA


40 bead populations and counting



- Map of 40 bead populations acquired on Mabtech EYRA
- Add assay-specific gates in Opal, linking bead populations to the associated analyte

Mabtech EYRA will save you time


Get your plates read before others even begin!





MABTECH EYRA - A SINGLE-PURPOSE CONFOCAL MICROSCOPE

- Own mechanics: X/Y and Z control for microscope motors
- Own electronics optimized for ultra-fast data-acquisition chain
- Laser-scanning technology
- Own high-speed FPGA electronics
- Patent-pending μm -accuracy autofocus technology



RAWsphere™ - MATHEMATICAL MODELING YIELDS ROBUST DATA

- Patent-pending RAWsphere™ algorithm
- Mabtech beads for ultra-precise data
- Mathematical modeling of bead fluorescence
- Inverse problems: Mathematical data extraction





Image data



Hypothesis after 2 iterations



Image data



Hypothesis after 4 iterations



Image data



Hypothesis after 8 iterations



Image data



Hypothesis after 16 iterations



Image data



Hypothesis after 32 iterations



Image data



Hypothesis after 64 iterations



Image data





Hypothesis after 100 iterations



RAWsphere™ - MATHEMATICAL MODELING YIELDS ROBUST DATA

- Inspired by the successful RAWspot™ algorithm for ELISpot and FluroSpot
- Resolves bead clusters effortlessly
- Impervious to bright-to-dim bead interactions
- Conceptualized in collaboration with Prof. Joakim Jaldén





Benefits of EYRAplex kits

	EYRAplex	Competitor
Magnetic bead-based assay	✓	✓
Mabtech's trusted mAb pairs	✓	✗
Goodbye heterophilic interference	✓	✗
Compatible with EYRA and most flow cytometers	✓	✗



Benefits of Mabtech EYRA


	EYRA	Competitor
Detect >30 analytes	✓	✓
Fluidics-free: no sheath fluid, no waste	✓	✗
Start, read and shutdown in <20 minutes	✓	✗
Export multiple plates into one Excel file	✓	✗
Worry-free coffee during readout	✓	✗

Assays & instruments designed for you

Your trusted source for accurate and sensitive immunoassays

ELISpot


IFN- γ



Each spot is the footprint
of a responding cell

FluoroSpot


IL-22, IFN- γ , IL-5, IL-17A



Multiplexed
ELISpot

FociSpot


Influenza foci



Virus-
infected cells

ELISA

IFN- γ



Analyte
concentration

EYRAplex

30+ analytes



Bead-based
multiplex

**IRIS 2
ASTOR 2**



EYRA



Mabtech EYRA

Multiplex made effortless

Q&A

